

# AIR SCIENCE FILTER REFERENCE SHEET



			Purair® Model P5-24	Purair® Model P5-36	Purair® Model P5-48	P10XL-CART-RC-P Air Science EDU-ADA Fume Hood	EDU-P15-ENCB Air Science EDU-Classic Fume Hood	EDU-M-40 Mobile Fume Hood on Cart	EDU-M-48 Mobile Fume Hood on Cart	EDU-M-60 Mobile Fume Hood on Cart	FLOW-24 Purair® FLOW™ Cabinet	FLOW-36 Purair® FLOW™ Cabinet	FLOW-48 Purair® FLOW™ Cabinet
ASTS-015	ACI Plus! 11-lb. Carbon Filter	Inorganic Acids	1 Required	1 Required	2 Required								
ASTS-030	HEPA Filter	Powders, Particulates	1 Required	1 Required	2 Required						1 Required		2 Required
ASTS-030U	ULPA filter FLOW-24 FLOW-48										1 Required		2 Required
ASTS-001	GP Plus! 11-lb. Carbon Filter	Organics, Solvents, Alcohols	1 Required	1 Required	2 Required								
ASTS-012	AMM 11 lb. filter	Amines, Ammonia	1 Required	1 Required	2 Required								
ASTS-013	FOR 11 lb. filter	Formaldehyde, Glutaraldehyde	1 Required	1 Required	2 Required								
ASTS-016	SUL 11 lb. filter	High MW mercaptans, Sulphur compounds	1 Required	1 Required	2 Required								
ASTM-001	GP Plus 22lb main filter	Alcohols, Solvents				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-002	ACR 22lb main filter	Radioactive Iodine				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-003	ACM 22lb main filter	Mercury				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-012	AMM 22lb main filter	Amines, Ammonia				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-013	FOR 22lb main filter	Formaldehyde, Glutaraldehyde				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-015	ACI Plus 22lb main filter	Inorganic Acids				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-016	SUL 22lb main filter	Sulfur compounds				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-017	CYN 22lb main filter	Hydrogen Cyanide				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-030	HEPA filter	Powders, Biologicals				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-100	EDU 22lb main filter	Acids, Amines, Organics				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-110	Multilayered 22lb main filter	Combination of two single layer filters				1 Required	1 Required	1 Required	1 Required	2 Required			
ASTM-200	GP/HEPA main filter	Organic vapors, powders				1 Required	1 Required	1 Required	1 Required	2 Required			
ASVLP536-030	HEPA filter FLOW-36											1 Required	
ASVLP536-030U	ULPA filter FLOW-36											1 Required	
ASTM-PRF	6 pack of Prefilters	Dust, Particulates				1 Required	1 Required	1 Required	1 Required	2 Required			
ASVLF-PRF	6 pack of Prefilters	Dust, Particulates									1 Required	1 Required	2 Required
ASTS-PRF (6 pack)	Pre-filter for particulates	Dust	1 Required	1 Required	2 Required								

