



The all-new Touch Screen temperature display and alarm module offers industry-leading product security real-time sample simulated temperatures are constantly displayed, logged and supported by audible and visual alarms, remote alarm contacts, a 4-20mA output and an access port. A USB port is provided for easy data transfer.



ABT-16S-TS



ABT-36S-TS

Standard features:
Keyed door locks.

Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
ABT-HC-12S-TS	12	Cycle	1 Swing Solid	3 Adj	25	29 3/4	63 3/4	R134a	3/8	6
ABT-16S-TS	16	Cycle	1 Swing Solid	3 Adj	25	29 3/4	77	R134a	3/8	6
ABT-23S-TS	23	Cycle	1 Swing Solid	4 Adj	26 7/8	34 3/4	79 3/4	R134a	3/8	6
ABT-HC-26S-TS	26	Cycle	1 Swing Solid	4 Adj	28 3/8	36 1/2	79 3/4	R134a	3/8	6
ABT-36S-TS	36	Cycle	2 Swing Solid	8 Adj	39 5/8	34 3/4	79 3/4	R134a	1/2	9
ABT-49S-TS	49	Cycle	2 Swing Solid	8 Adj	54	34 3/4	79 3/4	R134a	1/2	7
ABT-HC-72S-TS	72	Cycle	3 Swing Solid	12 Adj	81	34	79 3/4	R134a	3/4	10

- Swivel casters (adds 4" to height)
- Door ajar alarm
- Low battery alarm
- Operation indicator
- Alarm validation (audible, visual, and remote alarm)
- Alarm mute with ring back
- Remote alarm contacts
- Keyed door lock
- Magnetic door gasket for positive seal
- Heavy duty door handle

- LED interior lights are safety shielded and switch controlled
- Casters with footed brake
- Probe access port (3/8")
- White powder-coated steel exterior and interior
- CFC-free refrigerant (R134a), foam insulation and packaging
- Three year parts and labor warranty, plus an additional two year compressor parts warranty
- Quality certified with 100% unit temperature testing
- Cycle defrost
- UL/C-UL listed