

Premier Stainless Steel Refrigerators

Premier Stainless Steel Refrigerators

The Premier Stainless Steel Refrigerators feature a microprocessor temperature controller, allowing for precise temperature control, recovery and maintenance versus mechanical style thermostats.

Critical samples and supplies are security backed with an extensive package of temperature security and alarm features.



The Premier SS Refrigerators take a large step forward in performance from past solid door general purpose units; offering a Digital temperature display and alarm module coupled with a microprocessor temperature controller.



Stainless steel interior floor; stainless steel exterior.

Quality certified with 100% unit temperature testing.

ABT-SSP-49

Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
ABT-HC-SSP-23	23	Cycle	1 Solid	3 Adj	26 3/4	33	77 1/4	R134a	3/8	6
ABT-HC-SSP-49	49	Cycle	2 Solid	6 Adj	54	32 7/8	77 1/4	R134a	1/2	9

- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- Adjustable operating temperature range: 1°C – 10°C
- Keyed door lock
- Magnetic door gasket for positive seal
- Integrated door handle
- Swivel casters (adds 6 1/2" to height)
- Interior light

- Probe access port (3/8")
- Stainless steel interior floor
- Stainless steel exterior
- High density urethane foam cabinet insulation
- CFC-free refrigerant (R134a)
- Foam insulation and packaging
- Quality certified with 100% unit temperature testing
- Two year parts and labor warranty, plus an additional three year compressor parts warranty
- Cycle defrost
- UL/C-UL listed