



Key Features

- Microprocessor temperature controller
- Adjustable operating temperature range:
Refrigerator: 36°F – 46°F (2°C- 8°C)
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- Keyed door lock
- Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage
- Magnetic door gasket for positive seal
- Heavy duty door handle
- LED interior lights are safety shielded and switch controlled
- Casters with footed brake
- Probe access port (3/8")
- White powder-coated steel exterior and interior
- CFC-free refrigerant (R134a)
- Pharmacy refrigerator/freezer toolkit and temperature logs
- Quality certified with 100% unit temperature testing
- Cycle defrost
- Two year parts and labor warranty, plus an additional three year compressor parts warranty
- UL/C-UL listed

Standard Pharmacy Glass & Solid Door Refrigerators

The all-new Pharmacy **Standard Series Glass & Solid Door Pharmacy Refrigerators** are purpose built for pharmacy and medical applications. The new Standard Series Glass & Solid Door Pharmacy Refrigerators feature a digital microprocessor temperature controller allowing for the precise temperature management necessary for critical samples and supplies.

A forced air directional refrigeration system provides superior temperature uniformity and recovery after door openings. An automatic cycle-defrost design prevents product damaging temperature spikes.

The all-new temperature display and alarm module offers product security with a monitoring and alarm system. Real-time sample simulated temperatures are constantly displayed and supported by audible/visual high and low temperature alarms. Remote alarm dry contacts and an access port are provided to simplify connections to existing facility monitoring systems.



Pharmacy Standard Glass & Solid Door Refrigerators



Standard features include a microprocessor temperature controller and display for superior temperature control, verification, and recovery.

Now Featuring!

Refrigerator freestanding thermometer with 3 year certificate of calibration

C/F convertible temperature display.

Casters with footed brake



PH-ABT-S49G

Standard Glass Door Refrigerator

Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
PH-ABT-S12G	12	Cycle	1 Swing Glass	3 Adj	25	29 3/4	63 3/4	R134a	3/8	6
PH-ABT-S16G	16	Cycle	1 Swing Glass	3 Adj	25	29 3/4	77	R134a	3/8	6
PH-ABT-HC-S23G	23	Cycle	1 Swing Glass	4 Adj	26 7/8	34 3/4	79 3/4	R134a	3/8	6
PH-ABT-S26G	26	Cycle	1 Swing Glass	4 Adj	28 3/8	36 1/2	79 3/4	R134a	3/8	6
PH-ABT-S33G	33	Cycle	2 Slide Glass	8 Adj	39 5/8	33 1/4	79 3/4	R134a	3/8	6
PH-ABT-S36G	36	Cycle	2 Swing Glass	8 Adj	39 5/8	34 3/4	79 3/4	R134a	1/2	9
PH-ABT-HC-S47G	47	Cycle	2 Slide Glass	8 Adj	54	33 1/4	79 3/4	R134a	1/2	7
PH-ABT-HC-S49G	49	Cycle	2 Swing Glass	8 Adj	54	34 3/4	79 3/4	R134a	1/2	7
PH-ABT-HC-S69G	69	Cycle	3 Slide Glass	12 Adj	81	33 1/4	79 3/4	R134a	3/4	10
PH-ABT-HC-S72G	72	Cycle	3 Swing Glass	12 Adj	81	34	79 3/4	R134a	3/4	10

Standard Solid Door Refrigerator

Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
PH-ABT-HC-S12S	12	Cycle	1 Swing Solid	3 Adj	25	29 3/4	63 3/4	R134a	3/8	6
PH-ABT-S16S	16	Cycle	1 Swing Solid	3 Adj	25	29 3/4	77	R134a	3/8	6
PH-ABT-S23S	23	Cycle	1 Swing Solid	4 Adj	26 7/8	34 3/4	79 3/4	R134a	3/8	6
PH-ABT-HC-S26S	26	Cycle	1 Swing Solid	4 Adj	28 3/8	36 1/2	79 3/4	R134a	3/8	6
PH-ABT-S36S	36	Cycle	2 Swing Solid	8 Adj	39 5/8	34 3/4	79 3/4	R134a	1/2	9
PH-ABT-S49S	49	Cycle	2 Swing Solid	8 Adj	54	34 3/4	79 3/4	R134a	1/2	7
PH-ABT-HC-S72S	72	Cycle	3 Swing Solid	12 Adj	81	34	79 3/4	R134a	3/4	10