

Product Name:	ABT-HC-SLS-12
Description:	<p>Standard Solid Door Laboratory Refrigerator 12 Cu. Ft. capacity One swing solid door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door lock Heavy duty door handle Magnetic door gasket for positive seal Three adjustable shelves (adjustable in ½" increments) Forced draft air circulation LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (3/8") White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 25" W x 29 ¾" D x 65 ¾" H with casters; 62 ¼" H without casters Cycle defrost 115V, 60 Hz, 3 Amps, 1/5 HP ETL listed Shipping weight: 227 lbs.</p>

Product Name:	ABT-SLS-16
Description:	<p>Standard Solid Door Laboratory Refrigerator 16 Cu. Ft. capacity One swing solid door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Three adjustable shelves (adjustable in ½" increments) Heavy duty door handle Keyed door lock Magnetic door gasket for positive seal LED interior lights are safety shielded and switch controlled Casters with footed brake Probe access port (3/8") White powder-coated steel exterior and interior CFC-free refrigerant (R134a), foam insulation and packaging One year parts and labor warranty, plus an additional four year compressor parts warranty Quality certified with 100% unit temperature testing Exterior dimensions: 25"W x 29 ¾"D x 77"H Cycle defrost 115V, 60 Hz, 6 Amps, 3/8 HP UL/C-UL listed Shipping weight: 277 lbs.</p>

Product Name:	ABT-SLS-23
Description:	<p>Standard Solid Door Laboratory Refrigerator 23 Cu. Ft. capacity One swing solid door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle Four adjustable shelves (adjustable in ½" increments) LED interior lights are safety shielded and switch controlled Casters with footed brake Probe access port (3/8") White powder-coated steel exterior and interior CFC-free refrigerant (R134a), foam insulation and packaging One year parts and labor warranty, plus an additional four year compressor parts warranty Quality certified with 100% unit temperature testing Exterior dimensions: 26 7/8"W x 34 3/4"D x 79 3/4"H Cycle defrost 115V, 60 Hz, 6 Amps, 3/8 HP UL/C-UL listed Shipping weight: 321 lbs.</p>

Product Name:	ABT-HC-SLS-26
Description:	<p>Standard Solid Door Laboratory Refrigerator 26 Cu. Ft. capacity One swing solid door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: 1°C – 10°C Four adjustable shelves (adjustable in ½" increments) Forced draft air circulation Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (¾") White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 28 ¾" W x 36 ½" D x 81 ¾" H with casters; 78 ¼" H without casters Cycle defrost 115V, 60 Hz, 3 Amps, 1/5 HP ETL listed Shipping weight: 275 lbs.</p>

Product Name:	ABT-SLS-36
Description:	<p>Standard Solid Door Laboratory Refrigerator 36 Cu. Ft. capacity Two swing solid doors, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door locks Magnetic door gasket for positive seal Eight adjustable shelves (adjustable in ½" increments) Heavy duty door handles LED interior lights are safety shielded and switch controlled Casters with footed brake Probe access port (3/8") White powder-coated steel exterior and interior CFC-free refrigerant (R134a), foam insulation and packaging One year parts and labor warranty, plus an additional four year compressor parts warranty Quality certified with 100% unit temperature testing Exterior dimensions: 39 5/8"W x 34 3/4"D x 79 3/4"H Cycle defrost 115V, 60 Hz, 9 Amps, 1/2 HP UL/C-UL listed Shipping weight: 422 lbs.</p>

Product Name:	ABT-SLS-49
Description:	<p>Standard Solid Door Laboratory Refrigerator 49 Cu. Ft. capacity Two swing solid doors, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door locks Magnetic door gasket for positive seal Eight adjustable shelves (adjustable in ½" increments) Heavy duty door handles LED interior lights are safety shielded and switch controlled Casters with footed brake probe access port (3/8") White powder-coated steel exterior and interior CFC-free refrigerant (R134a), foam insulation and packaging One year parts and labor warranty, plus an additional four year compressor parts warranty Quality certified with 100% unit temperature testing Exterior dimensions: 54"W x 34 ¾"D x 79 ¾"H Cycle defrost 115V, 60 Hz, 7 Amps, 1/2 HP UL/C-UL listed Shipping weight: 528 lbs.</p>

Product Name:	ABT-HC-SLS-72
Description:	<p>Standard Solid Door Laboratory Refrigerator 72 Cu. Ft. capacity Three swing solid doors, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door locks Heavy duty door handles Magnetic door gasket for positive seal Twelve adjustable shelves (adjustable in ½" increments) Forced draft air circulation LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (¾") White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 81" W x 34" D x 81 ¾" H with casters; 78 ¼" H without casters Cycle defrost 115V, 60 Hz, 7.5 Amps, 1/2 HP ETL listed Shipping weight: 600 lbs.</p>