

Product Name:	ABT-HC-SSP-49
Description:	<p>Premier Stainless Steel Laboratory Refrigerator  49 Cu. Ft. Capacity  Two stainless steel swing doors  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range: 1°C -10°C  Six adjustable shelves  Forced draft air circulation  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Swivel casters (adds 6 ½" to height)  Interior light  Probe access port (3/8")  Stainless steel interior floor  Stainless steel exterior  High density urethane foam cabinet insulation  Hydrocarbon, natural refrigerant (R290)  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 54"W x 35"D x 84"H with casters; 77 ½"H without casters  Cycle defrost  115V, 60 Hz, 4.5 Amps, 1/4 HP  ETL listed  Shipping weight: 450 lbs.</p>