

Key Features

- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- Adjustable operating temperature range:
Refrigerator: 36°F – 46°F (2°C -8°C)
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable
- Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall
- Bottom footed shelf to allow for proper air circulation
- Back wall evaporator cover
- Forced draft circulation
- Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- LED interior light
- Probe access port (3/8")
- Powder coat paint exterior
- High density urethane foam cabinet insulation
- Hydrocarbon natural refrigerant (R600a)
- Pharmacy refrigerator/freezer toolkit and temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- Quality certified with 100% unit temperature testing
- Note: This unit must have 4" clearance on sides and back for adequate ventilation
- Cycle defrost
- Two year parts and labor warranty, plus an additional three year compressor parts warranty
- UL/C-UL listed

Premier Compact Pharmacy/Vaccine Refrigerators

An economical choice for low to mid volume pharmacies



Microprocessor temperature controller; Digital temperature display; Audible and visual high/low level alarms

Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall

PH-ABT-HC-10PG


Requires at least 4" of clearance on each side for adequate ventilation

Pharmacy Compact Refrigerators



Standard features include a microprocessor temperature controller and display for superior temperature control, verification, and rapid recovery after door openings.

Keyed door lock standard



Available in either solid or glass door

PH-ABT-HC-10PS

Now Featuring!
Refrigerator freestanding thermometer with 3 year certificate of calibration



Premier Refrigerator

Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
PH-ABT-HC-10PG	10.5	Cycle	1 Solid	6 Adj	23 3/4	26 1/2	59	R600a	1/3	2.3
PH-ABT-HC-10PS	10.5	Cycle	1 Glass	6 Adj	23 3/4	26 1/2	59	R600a	1/3	2.3