

Product Name:	ABT-HC-19-TS
Description:	<p>TempLog Premier Glass Door Laboratory Refrigerator</p> <p>19 Cu. Ft. capacity</p> <p>One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible</p> <p>Microprocessor temperature controller</p> <p>Adjustable operating temperature range: 1°C – 10°C</p> <p>Touch Screen Digital Temperature Display and Alarm Module</p> <p>Password-protected control panel</p> <p>Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</p> <p>Active temperature scale bar display</p> <p>Resettable 1-point min/max history</p> <p>Digital 24-hour temperature chart</p> <p>Temperature data logging</p> <p>USB port for data transfer and software updates</p> <p>Real-time operational display clock</p> <p>4-20mA output</p> <p>Electronic troubleshooting guide</p> <p>Electronic QR Code to access Owner's Manual & technical support contact Information</p> <p>Adjustable alarm volume</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Door ajar alarm</p> <p>Low battery alarm</p> <p>Operation indicator</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Remote alarm contacts</p> <p>Four adjustable shelves (adjustable in ½" increments)</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handle</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Swivel casters (adds 4" to height)</p> <p>Probe access port (3/8")</p> <p>White laminated vinyl exterior</p> <p>White aluminum liner and stainless steel floor interior</p> <p>HFC-free refrigerant (R290a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 27"W x 26 3/8"D x 82 5/8"H with casters; 78 5/8" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 5.4 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 345 lbs.</p>

Product Name:	ABT-HC-19-TS-LH
Description:	<p>TempLog Premier Glass Door Laboratory Refrigerator</p> <p>19 Cu. Ft. capacity</p> <p>One swing glass door, self-closing (energy efficient, double pane vacuum insulated), LEFT hinged, not reversible</p> <p>Microprocessor temperature controller</p> <p>Adjustable operating temperature range: 1°C - 10°C</p> <p>Touch Screen Digital Temperature Display and Alarm Module</p> <p>Password-protected control panel</p> <p>Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</p> <p>Active temperature scale bar display</p> <p>Resettable 1-point min/max history</p> <p>Digital 24-hour temperature chart</p> <p>Temperature data logging</p> <p>USB port for data transfer and software updates</p> <p>Real-time operational display clock</p> <p>4-20mA output</p> <p>Electronic troubleshooting guide</p> <p>Electronic QR Code to access Owner's Manual & technical support contact Information</p> <p>Adjustable alarm volume</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Door ajar alarm</p> <p>Low battery alarm</p> <p>Operation indicator</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Remote alarm contacts</p> <p>Four adjustable shelves (adjustable in ½" increments)</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handle</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Swivel casters (adds 4" to height)</p> <p>Probe access port (3/8")</p> <p>White laminated vinyl exterior</p> <p>White aluminum liner and stainless steel floor interior</p> <p>HFC-free refrigerant (R290a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 27"W x 26 3/8"D x 82 5/8"H with casters; 78 5/8" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 5.4 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 345 lbs.</p>

Product Name:	ABT-HC-23-TS
Description:	<p>TempLog Premier Glass Door Laboratory Refrigerator</p> <p>23 Cu. Ft. capacity</p> <p>One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible</p> <p>Microprocessor temperature controller</p> <p>Adjustable operating temperature range: 1°C - 10°C</p> <p>Touch Screen Digital Temperature Display and Alarm Module</p> <p>Password-protected control panel</p> <p>Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</p> <p>Active temperature scale bar display</p> <p>Resettable 1-point min/max history</p> <p>Digital 24-hour temperature chart</p> <p>Temperature data logging</p> <p>USB port for data transfer and software updates</p> <p>Real-time operational display clock</p> <p>4-20mA output</p> <p>Electronic troubleshooting guide</p> <p>Electronic QR Code to access Owner's Manual & technical support contact Information</p> <p>Adjustable alarm volume</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Door ajar alarm</p> <p>Low battery alarm</p> <p>Operation indicator</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Remote alarm contacts</p> <p>Four adjustable shelves (adjustable in ½" increments)</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handle</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Swivel casters (adds 4" to height)</p> <p>Probe access port (3/8")</p> <p>White laminated vinyl exterior</p> <p>White aluminum liner and stainless steel floor interior</p> <p>HFC-free refrigerant (R290a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 27"W x 31 3/8"D x 82 5/8"H with casters; 78 5/8" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 5.4 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 340 lbs.</p>

Product Name:	ABT-HC-26-TS
Description:	<p>TempLog Premier Glass Door Laboratory Refrigerator</p> <p>26 Cu. Ft. capacity</p> <p>One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible</p> <p>Microprocessor temperature controller</p> <p>Adjustable operating temperature range: 1°C - 10°C</p> <p>Touch Screen Digital Temperature Display and Alarm Module</p> <p>Password-protected control panel</p> <p>Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</p> <p>Active temperature scale bar display</p> <p>Resettable 1-point min/max history</p> <p>Digital 24-hour temperature chart</p> <p>Temperature data logging</p> <p>USB port for data transfer and software updates</p> <p>Real-time operational display clock</p> <p>4-20mA output</p> <p>Electronic troubleshooting guide</p> <p>Electronic QR Code to access Owner's Manual & technical support contact Information</p> <p>Adjustable alarm volume</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Door ajar alarm</p> <p>Low battery alarm</p> <p>Operation indicator</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Remote alarm contacts</p> <p>Four adjustable shelves (adjustable in ½" increments)</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handle</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Swivel casters (adds 4" to height)</p> <p>Probe access port (3/8")</p> <p>White laminated vinyl exterior</p> <p>White aluminum liner and stainless steel floor interior</p> <p>HFC-free refrigerant (R290a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 30"W x 31 3/8"D x 82 5/8" H with casters; 78 5/8" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 5.4 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 355 lbs.</p>

Product Name:	ABT-HC-47-TS
Description:	<p>TempLog Premier Glass Door Laboratory Refrigerator</p> <p>47 Cu. Ft. capacity</p> <p>Two sliding glass doors, self-closing (energy efficient, double pane vacuum insulated)</p> <p>Microprocessor temperature controller</p> <p>Adjustable operating temperature range: 1°C - 10°C</p> <p>Touch Screen Digital Temperature Display and Alarm Module</p> <p>Password-protected control panel</p> <p>Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</p> <p>Active temperature scale bar display</p> <p>Resettable 1-point min/max history</p> <p>Digital 24-hour temperature chart</p> <p>Temperature data logging</p> <p>USB port for data transfer and software updates</p> <p>Real-time operational display clock</p> <p>4-20mA output</p> <p>Electronic troubleshooting guide</p> <p>Electronic QR Code to access Owner's Manual & technical support contact Information</p> <p>Adjustable alarm volume</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Door ajar alarm</p> <p>Low battery alarm</p> <p>Operation indicator</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Remote alarm contacts</p> <p>Eight adjustable shelves (adjustable in ½" increments)</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Integrated door handles</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Swivel casters (adds 4" to height)</p> <p>Probe access port (3/8")</p> <p>White laminated vinyl exterior</p> <p>White aluminum liner and stainless steel floor interior</p> <p>HFC-free refrigerant (R290a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 54 1/8"W x 33 3/8"D x 82 5/8" H with casters; 78 5/8" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 9.6 Amps, 1/2 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 490 lbs.</p>

Product Name:	ABT-HC-49-TS
Description:	<p>TempLog Premier Glass Door Laboratory Refrigerator</p> <p>49 Cu. Ft. capacity</p> <p>Two swing glass doors, self-closing (energy efficient, double pane vacuum insulated)</p> <p>Microprocessor temperature controller</p> <p>Adjustable operating temperature range: 1°C - 10°C</p> <p>Touch Screen Digital Temperature Display and Alarm Module</p> <p>Password-protected control panel</p> <p>Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</p> <p>Active temperature scale bar display</p> <p>Resettable 1-point min/max history</p> <p>Digital 24-hour temperature chart</p> <p>Temperature data logging</p> <p>USB port for data transfer and software updates</p> <p>Real-time operational display clock</p> <p>4-20mA output</p> <p>Electronic troubleshooting guide</p> <p>Electronic QR Code to access Owner's Manual & technical support contact Information</p> <p>Adjustable alarm volume</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Door ajar alarm</p> <p>Low battery alarm</p> <p>Operation indicator</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Remote alarm contacts</p> <p>Eight adjustable shelves (adjustable in ½" increments)</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Heavy duty door handles</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Swivel casters (adds 4" to height)</p> <p>Probe access port (3/8")</p> <p>White laminated vinyl exterior</p> <p>White aluminum liner and stainless steel floor interior</p> <p>HFC-free refrigerant (R290a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 54 1/8"W x 31 3/8"D x 82 5/8" H with casters; 78 5/8" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 8.5 Amps, 1/2 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 470 lbs.</p>

Product Name:	ABT-HC-69-TS
Description:	<p>TempLog Premier Glass Door Laboratory Refrigerator</p> <p>69 Cu. Ft. capacity</p> <p>Three sliding glass doors, self-closing (energy efficient, double pane vacuum insulated)</p> <p>Microprocessor temperature controller</p> <p>Adjustable operating temperature range: 1°C - 10°C</p> <p>Touch Screen Digital Temperature Display and Alarm Module</p> <p>Password-protected control panel</p> <p>Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</p> <p>Active temperature scale bar display</p> <p>Resettable 1-point min/max history</p> <p>Digital 24-hour temperature chart</p> <p>Temperature data logging</p> <p>USB port for data transfer and software updates</p> <p>Real-time operational display clock</p> <p>4-20mA output</p> <p>Electronic troubleshooting guide</p> <p>Electronic QR Code to access Owner's Manual & technical support contact Information</p> <p>Adjustable alarm volume</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Door ajar alarm</p> <p>Low battery alarm</p> <p>Operation indicator</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Remote alarm contacts</p> <p>Twelve adjustable shelves (adjustable in ½" increments)</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Integrated door handles</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Swivel casters (adds 4" to height)</p> <p>Probe access port (3/8")</p> <p>White laminated vinyl exterior</p> <p>White aluminum liner and stainless steel floor interior</p> <p>HFC-free refrigerant (R290a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 78 1/8"W x 33 3/8"D x 82 5/8" H with casters; 78 5/8" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 6.5 Amps, 3/4 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 682 lbs.</p>

Product Name:	ABT-HC-72-TS
Description:	<p>TempLog Premier Glass Door Laboratory Refrigerator</p> <p>72 Cu. Ft. capacity</p> <p>Three swing glass doors, self-closing (energy efficient, double pane vacuum insulated)</p> <p>Microprocessor temperature controller</p> <p>Adjustable operating temperature range: 1°C - 10°C</p> <p>Touch Screen Digital Temperature Display and Alarm Module</p> <p>Password-protected control panel</p> <p>Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</p> <p>Active temperature scale bar display</p> <p>Resettable 1-point min/max history</p> <p>Digital 24-hour temperature chart</p> <p>Temperature data logging</p> <p>USB port for data transfer and software updates</p> <p>Real-time operational display clock</p> <p>4-20mA output</p> <p>Electronic troubleshooting guide</p> <p>Electronic QR Code to access Owner's Manual & technical support contact Information</p> <p>Adjustable alarm volume</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Door ajar alarm</p> <p>Low battery alarm</p> <p>Operation indicator</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Remote alarm contacts</p> <p>Twelve adjustable shelves (adjustable in ½" increments)</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Heavy duty door handles</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Swivel casters (adds 4" to height)</p> <p>Probe access port (3/8")</p> <p>White laminated vinyl exterior</p> <p>White aluminum liner and stainless steel floor interior</p> <p>HFC-free refrigerant (R290a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 78 1/8"W x 33 1/8"D x 82 5/8" H with casters; 78 5/8" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 8.5 Amps, 1/2 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 677 lbs.</p>