



Material Safety Data Sheet Thymolphthalein

Section 1 - Chemical Product and Company Identification

MSDS Name:

Thymolphthalein

Catalog Numbers:

LC26060

Synonyms:

None

Company Identification:

LabChem Inc

200 William Pitt Way

Pittsburgh, PA 15238

Company Phone Number:

(412) 826-5230

Emergency Phone Number:

(800) 424-9300

CHEMTREC Phone Number:

(800) 424-9300

Section 2 – Composition, Information on Ingredients

CAS#	Chemical Name:	Percent
125-20-2	Thymolphthalein	99

Section 3 - Hazards Identification

Emergency Overview

Appearance: White needles or crystals.**Danger!** May cause respiratory tract irritation. May cause skin irritation. May cause eye irritation. May cause digestive tract irritation. The toxicological properties of this material have not been fully investigated.**Target Organs:** None.

Potential Health Effects

Eye:

May cause eye irritation.

Skin:

May cause skin irritation.

Ingestion:

The toxicological properties of this substance have not been fully investigated.

Inhalation:

The toxicological properties of this substance have not been fully investigated.

Chronic:

Not available.



Material Safety Data Sheet Thymolphthalein

Section 4 - First Aid Measures

Eyes:

Flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Get medical aid.

Skin:

Get medical aid. Flush skin with plenty of water for at least 15 minutes while removing contaminated clothing and shoes. Wash clothing before reuse.

Ingestion:

Call a poison control center. If swallowed, do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Get medical aid.

Inhalation:

Remove from exposure and move to fresh air immediately. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical aid. Do NOT use mouth-to-mouth resuscitation.

Notes to Physician:

Treat symptomatically and supportively.

Section 5 - Fire Fighting Measures

General Information:

As in any fire, wear a self-contained breathing apparatus in pressure-demand, MSHA/NIOSH (approved or equivalent), and full protective gear. Dusts at sufficient concentrations can form explosive mixtures with air. During a fire, thermal decomposition or combustion may generate irritating and highly toxic gases.

Extinguishing Media:

Use water spray, dry chemical, or carbon dioxide, or chemical foam.

Autoignition Temperature:

Not available.

Flash Point:

Not available.

NFPA Rating:

(estimated) Health: 1; Flammability: 0; Instability: 0

Explosion Limits:

Not available.

Section 6 - Accidental Release Measures

General Information:

Use proper personal protective equipment as indicated in Section 8.

Spills/Leaks:

Sweep up or absorb material, then place into a suitable clean, dry, closed container for disposal. Avoid generating dusty conditions.



Material Safety Data Sheet Thymolphthalein

Section 7 - Handling and Storage

Handling:

Wash thoroughly after handling. Use with adequate ventilation. Minimize dust generation and accumulation. Avoid contact with skin and eyes. Keep container tightly closed. Avoid ingestion and inhalation.

Storage:

Keep container closed when not in use. Store in a tightly closed container. Store in a cool, dry, well-ventilated area away from incompatible substances.

Section 8 - Exposure Controls, Personal Protection

Engineering Controls:

Facilities storing or utilizing this material should be equipped with an eyewash facility and a safety shower. Use adequate ventilation to keep airborne concentrations low.

Exposure Limits:

Chemical Name:	ACGIH	NIOSH	OSHA
Thymolphthalein	none listed	none listed	none listed

OSHA Vacated PELs:

Thymolphthalein: No OSHA Vacated PELs are listed for this chemical.

Personal Protective Equipment

Eyes:

Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166.

Skin:

Wear appropriate gloves to prevent skin exposure.

Clothing:

Wear appropriate protective clothing to prevent skin exposure.

Respirators:

Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN 149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.

Section 9 - Physical and Chemical Properties

Physical State:	Solid
Color:	White needles or crystals
Odor:	None
pH:	Not available
Vapor Pressure:	Not available
Vapor Density:	Not available
Evaporation Rate:	Not available
Viscosity:	Not available
Boiling Point:	Not available
Freezing/Melting Point:	250°C
Decomposition Temperature:	Not available
Solubility in water:	Insoluble



Material Safety Data Sheet Thymolphthalein

Specific Gravity/Density: Not available
Molecular Formula: C₂₈H₃₀O₄
Molecular Weight: 94.2076

Section 10 - Stability and Reactivity

Chemical Stability:

Stable under normal temperatures and pressures.

Conditions to Avoid:

High temperatures, incompatible materials, and strong oxidizers.

Incompatibilities with Other Materials:

Strong oxidizers.

Hazardous Decomposition Products:

Carbon monoxide, carbon dioxide.

Hazardous Polymerization:

Has not been reported.

Section 11 - Toxicological Information

RTECS:

CAS# 125-20-2: unlisted.

LD50/LC50:

Not available.

Carcinogenicity:

CAS# 125-20-2: Not listed by ACGIH, IARC, NTP, or CA Prop 65.

Epidemiology:

No information found

Teratogenicity:

No information found

Reproductive:

No information found

Mutagenicity:

No information found

Neurotoxicity:

No information found

Section 12 - Ecological Information

No information found

Section 13 - Disposal Considerations

Dispose of in accordance with Federal, State, and local regulations.



Material Safety Data Sheet Thymolphthalein

Section 14 - Transport Information

US DOT

Shipping Name: Not regulated
Hazard Class:
UN Number:
Packing Group:

Section 15 - Regulatory Information

US Federal

TSCA:

CAS# 125-20-2 is listed on the TSCA inventory.

SARA Reportable Quantities (RQ):

None of the chemicals in this material have an RQ.

CERCLA/SARA Section 313:

No chemicals are reportable under Section 313.

OSHA - Highly Hazardous:

None of the chemicals in this product are considered highly hazardous by OSHA.

US State

State Right to Know:

CAS# 125-20-2 is not present on state lists from California, Pennsylvania, Minnesota, Massachusetts, Florida, or New Jersey.

California Regulations:

None

European/International Regulations

Canadian DSL/NDL:

CAS# 125-20-2 is listed on Canada's DSL List.

Canada Ingredient Disclosure List:

CAS# 125-20-2 is not listed on the Canadian Ingredient Disclosure List.

Section 16 - Other Information

MSDS Creation Date: November 15, 2007

Revision Date: December 7, 2007

Information in this MSDS is from available published sources and is believed to be accurate. No warranty, express or implied, is made and LabChem Inc. assumes no liability resulting from the use of this MSDS. The user must determine suitability of this information for his application.