

Product Name:	PH-ABT-HC-RFC7
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  7 Cu. Ft. capacity (5.2 Ref./1.5 Freezer)  Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  3 adjustable shelves refrigerator/2 adjustable shelves freezer, shelves have back rails to prevent product from being pushed against back wall  Footed bottom shelf refrigerator, and freezer  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  LED Interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 ¾ "W x 24 ¼ "D x 53 ½ "H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.3 Amps refrigerator, 1 Amps freezer, 1/6 HP refrigerator, 1/8 HP freezer  UL/C-UL listed  Shipping weight: 169 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC7A
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  7 Cu. Ft. capacity (5.2 Ref./1.7 Freezer)  Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  3 adjustable shelves refrigerator/2 adjustable shelves freezer, shelves have guard rails to prevent product from being pushed against back wall  Footed bottom shelf refrigerator, and freezer  Back wall evaporator cover  Forced draft circulation  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccines in storage Keyed door locks  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  LED Interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 3/4 "W x 24 1/4"D x 53 1/8 "H  Auto defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.3 Amps refrigerator, 1 Amps freezer, 1/6 HP refrigerator, 1/8 HP freezer  UL/C-UL listed  Shipping weight: 169 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC9
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  9 Cu. Ft. capacity (4.6 Ref./4.2 Freezer)  Two exterior solid swing doors right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall  Back wall evaporator cover, Refrigerator  Forced draft circulation Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handle, freezer; heavy duty door handle, refrigerator  Leveling legs  Stacking kit  LED interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 ¾ "W x 25 ½ "D x 66 ¾ "H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 235 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC9G
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  9 Cu. Ft. Capacity (4.6 Ref./4.2 Freezer)  Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  2 adjustable/1 fixed shelves refrigerator with guard rail on back to prevent product from being pushed against back wall/2 fixed shelves freezer  Back wall evaporator cover, Refrigerator  Forced draft circulation, Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccines in storage  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handle, freezer; heavy duty door handle, refrigerator  Leveling legs  Stacking kit  LED interior light, Refrigerator  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Exterior dimensions: 23 ¾" W x 25 ½"D x 66 ¾"H  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 268 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC9G-LH
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  9 Cu. Ft. Capacity (4.6 Ref./4.2 Freezer)  Two exterior swing doors (glass door refrigerator/solid door freezer), LEFT HINGED  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C - -25°C  2 adjustable/1 fixed shelves refrigerator with guard rail on back to prevent product from being pushed against back wall/2 fixed shelves freezer  Back wall evaporator cover, Refrigerator  Forced draft circulation, Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccines in storage  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handle, freezer; heavy duty door handle, refrigerator  Leveling legs  Stacking kit  LED interior light, Refrigerator  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Exterior dimensions: 23 ¾" W x 25 ½"D x 66 ¾"H  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 268 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC9-LH
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  9 Cu. Ft. capacity (4.6 Ref./4.2 Freezer)  Two exterior solid swing doors, left hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall  Back wall evaporator cover, Refrigerator  Forced draft circulation Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handle, freezer; heavy duty door handle, refrigerator  Leveling legs  Stacking kit  LED interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 ¾ "W x 25 ½ "D x 66 ¾ "H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 235 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC9SS
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  9 Cu. Ft. Capacity (4.6 Ref./4.2 Freezer)  Two stainless steel exterior swing doors, right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall  Back wall evaporator cover, Refrigerator  Forced draft circulation Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handles  Leveling legs  Stacking kit  LED interior lights  Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer)  Stainless steel exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 ¾" W x 24 ½"D x 66 ¾"H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 230 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC9SS-LH
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination 9 Cu. Ft. Capacity (4.6 Ref./4.2 Freezer)</p> <p>Two stainless steel exterior swing doors, left hinged, not reversible</p> <p>Microprocessor temperature controller, refrigerator and freezer</p> <p>Digital temperature display, refrigerator and freezer</p> <p>Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)</p> <p>Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable</p> <p>Audible and visual high and low temperature alarms, refrigerator and freezer</p> <p>Remote alarm contacts, refrigerator and freezer</p> <p>2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall</p> <p>Back wall evaporator cover, Refrigerator</p> <p>Forced draft circulation Refrigerator</p> <p>Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage</p> <p>Keyed door locks</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handles</p> <p>Leveling legs</p> <p>Stacking kit</p> <p>LED interior lights</p> <p>Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer)</p> <p>Stainless steel exterior</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Pharmacy refrigerator/freezer toolkit and temperature logs</p> <p>Quality certified with 100% unit temperature testing</p> <p>Two year parts and labor warranty, plus an additional three year compressor parts warranty</p> <p>Exterior dimensions: 23 ¾" W x 24 ½"D x 66 ¾"H</p> <p>Manual defrost freezer / cycle defrost refrigerator</p> <p>115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer</p> <p>UL/C-UL listed</p> <p>Shipping weight: 230 lbs.</p> <p>*Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC12
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  12 Cu. Ft. capacity (10.5 Ref./1.7 Freezer)  Two exterior solid swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Refrigerator  Two adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Freezer  Bottom footed shelf to allow air circulation, Refrigerator and Freezer  Forced draft circulation Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  LED interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 3/4" W x 26 1/2" D x 80 1/4" H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 2.3 Amps refrigerator, 1 Amp freezer, 1/3 HP refrigerator, 1/8 HP freezer  UL/C-UL listed  Shipping weight: 261 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC12A
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  12 Cu. Ft. capacity (10.5 Ref./1.3 Freezer)  Two exterior solid swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Refrigerator  Two adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Freezer  Bottom footed shelf to allow air circulation, Refrigerator and Freezer  Forced draft circulation Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccines in storage  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  LED interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 3/4" W x 26 1/2" D x 80 1/4" H  Auto defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 2.3 Amps refrigerator, 1 Amp freezer, 1/3 HP refrigerator, 1/8 HP freezer  UL/C-UL listed  Shipping weight: 261 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC12G
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  12 Cu. Ft. capacity (10.5 Ref./1.7 Freezer)  Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Refrigerator  Two adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Freezer  Bottom footed shelf to allow air circulation, Refrigerator and Freezer  Forced draft circulation Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  LED interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 3/4" W x 26 1/2" D x 80 1/4" H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 2.3 Amps refrigerator, 1 Amp freezer, 1/3 HP refrigerator, 1/8 HP freezer  UL/C-UL listed  Shipping weight: 261 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-HC-RFC12GA
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  12 Cu. Ft. capacity (10.5 Ref./1.3 Freezer)  Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Two Refrigerator/Freezer freestanding thermometers with 3 year certificate of calibration, with min/max memory, °F/°C switchable  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Refrigerator  Two adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Freezer  Bottom footed shelf to allow air circulation, Refrigerator and Freezer  Forced draft circulation Refrigerator  Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccines in storage  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  LED lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 ¼" W x 26 ½" D x 80 2/5" H  Auto defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 2.3 Amps refrigerator, 1 Amp freezer, 1/3 HP refrigerator, 1/8 HP freezer  UL/C-UL listed  Shipping weight: 261 lbs.  *Please note this unit requires 2 outlets.</p>

Product Name:	PH-ABT-RFC-16A
Description:	<p>Premier Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  16 cu. ft. capacity (13 cu. ft. refrigerator/3 cu. ft. freezer)  Two exterior solid swing, double insulated doors, right hinged, not reversible  Individual microprocessor temperature controllers for both the refrigerator and freezer  Independent temperature displays  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C - 8°C); Freezer: 5°F – 13°C (-15°C - -25°C)  Freestanding thermometer with 2 year certificate of calibration  Sample temperature in freezer does not rise above -20°C during defrost cycle  Both chambers feature forced draft circulation / fan driven air circulation  Uniformity of samples throughout the freezer (<math>\pm 2^{\circ}\text{C}</math> throughout chamber)  Uniformity of samples throughout the refrigerator (<math>\pm 1^{\circ}\text{C}</math> throughout chamber)  36 hour battery back-up on alarms, high/low temperature &amp; loss of power  Two isolated audible &amp; visual high/low alarm systems remote alarm contacts  Remote alarm contacts, refrigerator &amp; freezer  3 adjustable, vented stainless steel shelves refrigerator/ 1 fixed shelf freezer  Cord label stating Warning! Do not unplug the appliance or break circuit.  Expensive vaccines in storage  Leveling legs  Keyed door locks, refrigerator &amp; freezer  Magnetic door gasket for positive seal  Heavy duty door handles  LED interior lights are safety shielded  Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer)  White powder coated exterior  CFC-free refrigerant (R134a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 29 1/2"W x 34 3/4"D x 78"H  Auto defrost freezer / cycle defrost refrigerator  115 V, 60Hz, 7.5 Amps, 1/8 HP  UL/C-UL listed  Shipping weight: 291 lbs.</p>