

Product Name:	ABT-MFP-14
Description:	<p>Premier Manual Defrost Laboratory Freezer</p> <p>14 Cu. Ft. capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Four Inner "fast-freeze compartments" with doors for temperature maintenance and recovery</p> <p>Microprocessor Temperature Controller</p> <p>Adjustable operating temperature range: -15°-25°C</p> <p>Standard Digital Temperature Display and Alarm Module</p> <p>Sample simulated digital temperature display; .1° Resolution</p> <p>Adjustable temperature display in Celsius or Fahrenheit</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Resettable 1-point Min/Max History</p> <p>Remote alarm contacts</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Three fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Heavy duty door handle</p> <p>Caster platform (adds 4 1/2" to height)</p> <p>One probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R134a), foam insulation and packaging</p> <p>Two year parts and labor warranty; plus additional three year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: W 29 3/4" x 32 6/16 "D x 62 1/2" H without caster platform; 67" H with caster platform</p> <p>Manual defrost</p> <p>115V, 60 Hz, 5 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 189 lbs.</p>

Product Name:	ABT-MFP-17
Description:	<p>Premier Manual Defrost Laboratory Freezer</p> <p>17 Cu. Ft. capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Four Inner "fast-freeze compartments" with doors for temperature maintenance and recovery</p> <p>Microprocessor Temperature Controller</p> <p>Adjustable operating temperature range: -15°-25°C</p> <p>Standard Digital Temperature Display and Alarm Module</p> <p>Sample simulated digital temperature display; .1° Resolution</p> <p>Adjustable temperature display in Celsius or Fahrenheit</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Resettable 1-point Min/Max History</p> <p>Remote alarm contacts</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Four fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Heavy duty door handle</p> <p>Caster platform (adds 4 1/2" to height)</p> <p>One probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R134a), foam insulation and packaging</p> <p>Two year parts and labor warranty; plus additional three year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 34" W x 30 7/16" D x 67 3/4" H without caster platform; 72 1/4" with caster platform</p> <p>Manual defrost</p> <p>115V, 60 Hz, 5 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 203 lbs.</p>

Product Name:	ABT-MFP-20
Description:	<p>Premier Manual Defrost Laboratory Freezer</p> <p>20 Cu. Ft. capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Five inner "fast-freeze compartments" with doors for temperature maintenance and recovery</p> <p>Microprocessor Temperature Controller</p> <p>Adjustable operating temperature range: -15°-25°C</p> <p>Standard Digital Temperature Display and Alarm Module</p> <p>Sample simulated digital temperature display; .1° Resolution</p> <p>Adjustable temperature display in Celsius or Fahrenheit</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Resettable 1-point Min/Max History</p> <p>Remote alarm contacts</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Four fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Heavy duty door handle</p> <p>Caster platform (adds 4 1/2" to height)</p> <p>One probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R134a), foam insulation and packaging</p> <p>Two year parts and labor warranty; plus additional three year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 34" W x 32 9/16" D x 73 1/4 "H without caster platform; 77 3/4" H with caster platform</p> <p>Manual defrost</p> <p>115V, 60 Hz, 5 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 221 lbs.</p>

Product Name:	ABT-MFP-30
Description:	<p>Premier Manual Defrost Laboratory Freezer</p> <p>30 Cu. Ft. capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Five inner "fast-freeze compartments" with doors for temperature maintenance and recovery</p> <p>Microprocessor Temperature Controller</p> <p>Adjustable operating temperature range: -15°-25°C</p> <p>Standard Digital Temperature Display and Alarm Module</p> <p>Sample simulated digital temperature display; .1° Resolution</p> <p>Adjustable temperature display in Celsius or Fahrenheit</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Resettable 1-point Min/Max History</p> <p>Remote alarm contacts</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Four fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Stainless steel door handle</p> <p>Caster platform (adds 4 1/2" to height)</p> <p>One probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R404a), foam insulation and packaging</p> <p>Two year parts and labor warranty; plus additional three year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: W 35 3/4" x 34 1/4" D x 72 1/2" H without caster platform; 77" H with caster platform</p> <p>Manual defrost</p> <p>115V, 60 Hz, 6 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 380 lbs.</p>

Product Name:	ABT-MFP-2030
Description:	<p>Premier Manual Defrost Laboratory Freezer</p> <p>20 Cu. Ft. capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Five inner "fast-freeze compartments" with doors for temperature maintenance and recovery</p> <p>Microprocessor Temperature Controller</p> <p>Adjustable operating temperature range: -27°C – -33°C</p> <p>Standard Digital Temperature Display and Alarm Module</p> <p>Sample simulated digital temperature display; .1° Resolution</p> <p>Adjustable temperature display in Celsius or Fahrenheit</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Resettable 1-point Min/Max History</p> <p>Remote alarm contacts</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Four fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Heavy duty door handle</p> <p>Caster platform (adds 4 1/2" to height)</p> <p>One probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R134a), foam insulation and packaging</p> <p>Two year parts and labor warranty; plus additional three year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: 34" W x 32 9/16" D x 73 1/4" H without caster platform; 77 3/4" H with caster platform</p> <p>Manual defrost</p> <p>115V, 60 Hz, 5 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 221 lbs.</p>

Product Name:	ABT-MFP-3030
Description:	<p>Premier Manual Defrost Laboratory Freezer</p> <p>30 Cu. Ft. capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Five Inner "fast-freeze compartments" with doors for temperature maintenance and recovery</p> <p>Microprocessor Temperature Controller</p> <p>Adjustable operating temperature range: -27°-33°C</p> <p>Standard Digital Temperature Display and Alarm Module</p> <p>Sample simulated digital temperature display; .1° Resolution</p> <p>Adjustable temperature display in Celsius or Fahrenheit</p> <p>Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Resettable 1-point Min/Max History</p> <p>Remote alarm contacts</p> <p>Alarm validation (audible, visual, and remote alarm)</p> <p>Alarm mute with ring back</p> <p>Four fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Stainless steel door handle</p> <p>Caster platform (adds 4 1/2" to height)</p> <p>One probe access port (3/8")</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R404a), foam insulation and packaging</p> <p>Two year parts and labor warranty; plus additional three year compressor parts warranty</p> <p>Quality certified with 100% unit temperature testing</p> <p>Exterior dimensions: W 35 3/4 " 34 1/4" D x 72 1/2 "H without caster platform; 77" H with caster platform</p> <p>Manual defrost</p> <p>115V, 60 Hz, 6 Amps, 1HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 380 lbs.</p>