

Rocker 35 EZ

The Rocker 35 EZ is the next generation large capacity rocker from Labnet. This upgraded Rocker 35 is designed for durability in laboratories requiring a larger platform area and precise control over rocking parameters. It is backed by a 5 year warranty.

Both speed and time are easily selectable and the analog control has been upgraded to display discreet values. The tilt is now easily adjustable with knob on side. (See figure below). This system can easily handle up to 11 pounds or 5 kg of load so a double platform can be fully loaded without any variation to speed.

This unit comes in both single and dual platform models and includes a ribbed no-slip rubber mat to keep containers in place. This unit is also compatible with our Sticky Mats (I-5240) for applications requiring greater sample adhesion to the platform.



Adjustment knob

SPECIFICATIONS

| | |
|-------------------------|---|
| Speed range | 10 - 120 rpm |
| Timer | 10 - 120 min or continuous |
| Tilt range | Adjustable, $\pm 20^\circ$ |
| Maximum Load | 11 lb/5kg |
| Ambient operating range | -10° to 70°C , up to 85%RH |
| Dimensions (W x D x H) | |
| Single platform | 14.2 x 14.2 x 6.7 in / 36 x 36 x 17 cm |
| Double platform | 14.2 x 14.2 x 10.2 in / 36 x 36 x 26 cm |
| Weight | |
| Single platform | 14.3 lb/6.5.7 kg |
| Double platform | 17.2 lb/7.8kg |

| CAT NO. | DESCRIPTION |
|-----------|-------------------------------|
| S2035-E* | Rocker 35 EZ, single platform |
| S2035-DE* | Rocker 35 EZ, double platform |

*To order 230V units add -230V to the end of the catalog number. includes EU and UK cords.